

TECHNICAL DATA SHEET

002561-TDS-ENG-2026

CAPSULAS 3 INCOLORAS		
DESCRIPTION DCI: ---		DESCRIPTION DOE: ---
CAS Nº: ---	EC Nº: ---	AEMPS CODE: ---
MOL. WEIGHT: ---	MOL. FORMULA: ---	ARTICLE CODE: 002561

ATTRIBUTES	SHOULD BE
Chemical tests	
Disintegration test	< 15 min
Microbiological control	
TAMC	< 1000 CFU/g
Salmonella	Negative/10g
Escherichia coli	Negative/1g
Staphylococcus aureus	Negative/1g
Pseudomonas aeruginosa	Negative/1g
Dimensional and physical tests	
Diameter cap	5.86 - 5.84 mm
Body diameter	5.59 - 5.57 mm
Cap length	8.08 - 7.84 mm
Length body	13.74 - 13.45 mm
	129 - 114 µm
	150 - 109 µm
Weight	53.0 - 46.3 mg
Moisture	16.0 - 13.0 %
Capsule colour formulation	
Cap	Water (14.5%) Target Moisture Gelatin
Body	Water (14.5%) Target Moisture Gelatin

COMPLIES WITH

Manufacturer Specifications

STORAGE

Keep the containers tightly closed. Store in a fresh, cool and well-ventilated place.

REMARKS

The design weight of a "size 3" capsule is 50.0 mg, with the cap accounting for 40% and the body for 60% of the capsule's weight. The average capsule weight calculated from 100 capsules may vary within the range of 46.3-53.8 mg.

Raw material

GELATIN (**): Complies with the requirements of the current Editions of Eur. Ph. and USP/NF. It is of purely bovine origin, it complies with the current revision of the European guideline EMEA/410/01. Each supplier has a CEP for their product (R1-CEP-2000-027-Rev.02), R1-CEP-2000-029-Rev.05, R1-CEP-2000-045-Rev.03, R1-CEP- 2002-110-Rev.00 and R1-CEP-2001-211.Rev.01)

COLORANTS: They comply with Directive 2009/35, Commission Regulation 231/2012 and, where applicable, with the requirements of the Eur. Ph. and USP/NF pharmacopoeias.

TECHNICAL DATA SHEET

002561-TDS-ENG-2026

CAPSULAS 3 INCOLORAS		
DESCRIPTION DCI: ---		DESCRIPTION DOE: ---
CAS Nº: ---	EC Nº: ---	AEMPS CODE: ---
MOL. WEIGHT: ---	MOL. FORMULA: ---	ARTICLE CODE: 002561

PRINTING INKS: They comply with pharmaceutical regulations.

Capsules

It does not contain preservatives and has not been treated with ethylene oxide. They comply with the CPMP/ICH/283/95 guideline of the European Agency for the Evaluation of Medicinal Products (EMA) and with the European Pharmacopoeia for residual solvents.

DISINTEGRATION (***): Less than 15 minutes per test according to Eur. Ph.. They have a breakdown time of less than 5 minutes according to US Federal Specification 285a on the acid solubility test.

With this we certify that the above information has been approved by technical management in accordance with the specifications, as described by the applicable regulatory requirements.