

Ph. Eur. reference spectrum of mupirocin

Batch 1

1. Identification

Catalogue code: **M3805000**

2. Scientific Information

2.1 Intended use

Established for use with the monograph(s): **1450**

This official Reference Spectrum is to be used with the currently valid European Pharmacopoeia monograph(s) and/or general chapter(s) describing its suitable use.

Further information about the Reference Standards is available in the on-line catalogue currently at <http://crs.edqm.eu>.

2.2. Preparation

This spectrum has been recorded using 50 µl of a 20 g/l solution of mupirocin in methylene chloride placed on a disc of 300 mg of potassium bromide R after evaporation of the solvent. The absorption maxima of polystyrene at 2850 cm⁻¹ (3.51 µm), 1601 cm⁻¹ (6.25 µm) and 1028 cm⁻¹ (9.73 µm) have been superimposed on the spectrum.

3. Adoption


The suitability for intended use has been officially adopted by the European Pharmacopoeia Commission.

4. Signature

This document is approved by:

Ms Caroline Offerlé
Head of the Quality and Risk Management Section

EUROPEAN DIRECTORATE FOR THE QUALITY OF MEDICINES



copyright Council of Europe, 2001, all rights reserved

Mupirocin

ref: 11728