

CERTIFICATE OF ANALYSIS

PRODUCT ARGIRELINE® peptide solution C
CODE P06-PD010
BATCH F764/23

ANALYTICAL BULLETIN

TEST	SPECIFICATION	RESULT	METHOD
Physical appearance	Translucent solution	Correct	CC-10-01-B6
Colour	Colourless	Correct	CC-10-01-B6
Specific gravity	-	-	
pH	3.0-6.0	4.9	CC-10-01-B1
Viscosity (cPs) ¹	-	-	
Particle size (mm)	-	-	
Refractive index	-	-	
Dry residue (%)	-	-	
Identification (HPLC)	Similar retention time as in house ref.standard	Correct	CC-10-01-B19

MICROBIAL TEST

Total aerobic microbial count	<100 cfu/g	Correct	MI-10-01-D09
Total yeast and mould count	<10 cfu/g	Correct	MI-10-01-D09
<i>Escherichia coli</i>	Absence/g	Correct	MI-10-01-D09
<i>Pseudomonas aeruginosa</i>	Absence/g	Correct	MI-10-01-D09
<i>Staphylococcus aureus</i>	Absence/g	Correct	MI-10-01-D09
<i>Candida albicans</i>	Absence/g	Correct	MI-10-01-D09

Manufacturing date: 12/23
Retest date: 12/25

Released on: 19/12/2023


QUALITY CONTROL
Ana Poza

QUALITY ASSURANCE
Dr. Fernando Fernández

Printed on: 23/04/2024

FG-10-01/2 v.3

Sample Information
Sample Name : ARGIRELINE® peptide solution C
ID : P06-PD010 F764/23
Data Filename : 12122023P06-PD010.003.lcd
Date Acquired : 12/12/2023 11:12:19
Column : Kromasil 100-5C8, 250 x 4.6 mm



PeakTable

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Name
1	11.728	5705794	Argireline ®
Total			