

Qualicaps Europe (\$34.1). 28108 Alcobendas (Madrid) Telefono: +34 - 92 663 08 00 65x1434 - 97 (63 08 00)

## CERTIFICATE OF ANALYSIS

This is to certify that the information below has been approved by QA as conforming to Qualicaps Europe specifications, as described in the current technical manual and applicable regulatory requirements

**CUSTOMER:** 

**ITEM DESCRIPTION: 1AJAAJA** 

CUSTOMER ITEM: E01334F

LOT NUMBER:

E1209192 LOT QUANTITY APPROVED: 1080000

PRINTING: INK COLOR: NOT APPLICABLE

NOT APPLICABLE

COUNTRY: SPAIN

4900004513 CUSTOMER ORDER:

NAME:

APPROVAL DATE:

17/12/2012

**MANUFACTURING DATE: 10/12/2012** EXPIRY DATE:

**5 YEARS** 

GELATIN: Complies with the requirements of current EP and USP/NF editions.

It is of pure bovine origin, complies with revision in force of European Guideline EMEA/410/01. Each supplier has a Ph. Eur certificate of suitability for their product (R1-CEP-2000-027-Rev.02, R1-CEP-2000-029-Rev.04, R1-CEP-2002-110-Rev.00, R1-CEP-2001-211.Rev.01) and has confirmed that it is manufactured from raw material not imported from UK, Ireland and Portugal.

COLORANTS: Comply with the EC directives 2009/35 & 2008/128 and where applicable with the requirements of the European Pharmacopoeia and the FDA.

RINTING INKS: Comply with pharmaceutical regulations.

## FINISHED PRODUCT TESTS

CAPSULES: Do not contain preservatives and have not been treated with Ethylene Oxide. They comply with The European Agency for the Evaluation of Medicinal Products (EMEA) guideline CPMP/ICH/283/95 and the European Pharmacopoeia for residual solvents. DISINTEGRATION: Less than 15 minutes by the Ph. Eur. test and have a rupture time of less than 5 minutes by the USA Federal Specification 285A

PARAMETERS	LIMITS	REFERENCE	RESULT					
MICROBIOLOGICAL AND CHEMICAL TEST								
Total Aerobic Microbial Count (	TAMC)103 CELI/G	EP/USP	10 (Average)					
Salmonella	Negative / 10 G	EP/USP	Meets test					
E. Coli	Negative / 1 G	EP/USP	Meets test					
S. Aureus	Negative / 1 G	EP/USP	Meets test					
Ps. Aeruginosa	Negative / 1 G	EP/USP	Meets test					

DIAMETER	CAP (mm)	Meets spec	TM*	6,91	-	6,91	(max-min)
		Meets spec		6,61	-	6,61	(max-min)
ENGTH	CAP (mm)	Meets spec	TM*	9,9	-	9,7	(max-min)
SNOTH	BODY (mm)	Meets spec	TM*	16,9	-	16,5	(max-min)
END THICKNESS	CAP (micr)	Meets spec	QS*	176	-	113	(max-min)
BIND TIMOTENEOU	BODY (micr)	Meets spec		159	-	101	(max-min)
WEIGHT (MG)	2021 (mm-1,	Meets spec	TM*		(Average)		
MOISTURE (%)	***************************************	Meets spec		14,7			

<sup>\*</sup>TM = Qualicaps Technical Manual



<sup>\*</sup>QS = Qualicaps internal specification
"Meets test" means conforms to specification

CUSTOMER

LOT NUMBER: E1209192

**COUNTRY:** SPAIN

LOT QUANTITY APPROVED:

1080000

## CAPSULE COLOR FORMULATION SHEET

The design weight of a size 1 capsules is 76 mg, the cap representing 40% and the body 60% of the capsule weight. The mean capsule weight of 100 capsules can vary within the range of 70-82 mg.

CAP:

CODE

TITLE

**GM OF DYE PER 100 GM CAPSULE** 

AJA

**OPAQUE WHITE** 

TITANIUM DIOXIDE

2,000

Water (14,5%) target moisture

Gelatin

In the cap (40%)

30,4mg Corresponds to:

0,608

TITANIUM DIOXIDE

4,408 25,384 Water (14,5%) target moisture

Gelatin

BODY:

CODE

TITLE

**GM OF DYE PER 100 GM CAPSULE** 

AJA

**OPAQUE WHITE** 

TITANIUM DIOXIDE

2,000

Water (14,5%) target moisture

Gelatin

In the body (60%)

45,6mg Corresponds to:

0,912

TITANIUM DIOXIDE

6,612

Water (14,5%) target moisture

38,076

Gelatin

COLORANT TITANIUM DIOXIDE

C.I. NUMBER

E.E.C. NUMBER

## VISUAL AND PRINT QUALITY

All the capsules are controlled statistically to ensure conformance to specifications. These values are derived from the MIL.STD.105E. The quality specifications are as follow:

VISUAL DEFECTS AQL (%)

Major A 0,010 Major A Major B 0,040

Minor 0,25

PRINTING DEFECTS AQL (%)

0,010

Major B 0,040

Minor 1,0

QUALITY ASSURANCE MANAGEMENT **ENRIQUE REQUEJO** 

my toll

of 2

DATE:

17/12/2012