

Qualicaps Europe, S.A.U. Avda. Monte Valdelatas, 4 28108 Alcobendas, Madrid, Spain T +34 91 663 08 00 F +34 91 663 08 30



# QUALI-G™ GELATIN CAPSULE

## **CERTIFICATE OF ANALYSIS**

**CUSTOMER:** 

STOCK

ITEM DESCRIPTION: 1AJAAJA

**CUSTOMER ITEM: E66008F** 

LOT NUMBER:

E1703177

PRINTING: INK COLOR:

**LOT QUANTITY APPROVED: 3100000 NOT APPLICABLE NOT APPLICABLE**  COUNTRY:

France

SIZE:

**CUSTOMER ORDER:** 

NAME:

APPROVAL DATE:

EXPIRY DATE:

31/08/2017

MANUFACTURING DATE: 26/07/2017

5 YEARS

FINISHED PROD	OUCT TESTS
---------------	------------

PARAMETERS		LIMITS	REFERENCE	RESULT
V T - THE PART OF THE	er tract tale also go 2 to 12	MICROBIOLOGICAL AN	ID CHEMICAL TEST	
Total Aerobic Microbial Count (TAMC)		103 CFU/G	EP/USP	10 (Average)
Salmonella		Negative / 10 G	EP/USP	Meets test
E. Coli		Negative / 1 G	EP/USP	Meets test
S. Aureus		Negative / 1 G	EP/USP	Meets test
Ps. Aeruginosa		Negative / 1 G	EP/USP	Meets test
		DIMENSIONAL AND BUILD	(Olo )   Too	
Diameter	DIMENSIONAL AND PHYSICAL			0.04
Diameter	Cap (mm)	Meets spec Meets spec	TM* TM*	6,94 - 6,92 (max-i
ength	Body (mm) Cap (mm)	Meets spec	TM*	9,8 - 9,4 (max-i
Longui	Body (mm)	Meets spec	TM*	16,6 - 16,4 (max-r
End thickness	Cap (micr)	Meets spec	QS*	169 - 116 (max-r
	Body (micr)	Meets spec	QS*	185 - 118 (max-r
Neight (mg)		Meets spec	TM*	79 (Average)
Moisture (%)		Meets spec	TM*	14,4 (Average)
				×
****		VISUAL ATTRIBUTES	. 44.	
VISUAL DEFECTS AQL (%)		Major A Major B Minor 0,010 0,040 0,25	TM*	Meets
PRINTING DEFECTS AQL (%)		Major A Major B Minor 0,010 0,040 1,0	TM*	Meets (if applicable)

\*TM = Qualicaps Technical Manual

\*QS = Qualicaps Internal specification "Meets test" means conforms to specification



CUSTOMER:

STOCK

ITEM DESCRIPTION: 1AJAAJA

**CUSTOMER ITEM: E66008F** 

LOT NUMBER:

E1703177

**LOT QUANTITY APPROVED: 3100000** PRINTING:

**NOT APPLICABLE** 

INK COLOR:

**NOT APPLICABLE** 

COUNTRY:

France

SIZE:

**CUSTOMER ORDER:** 

NAME:

APPROVAL DATE:

31/08/2017

MANUFACTURING DATE: 26/07/2017

EXPIRY DATE:

5 YEARS

### CAPSULE COLOR FORMULATION

The design weight of a size 1 capsules is 76 mg, the cap representing 40% and the body 60% of the capsule weight. The mean capsule weight of 100 capsules can vary within the range of 70-82 mg.

CAP:

CODE

AJA

TITLE

GM OF DYE PER 100 GM CAPSULE

**OPAQUE WHITE** 

TITANIUM DIOXIDE

2.000

WATER (14,5%) TARGET MOISTURE

**GELATIN** 

BODY:

CODE

TITLE

GM OF DYE PER 100 GM CAPSULE

AJA

**OPAQUE WHITE** 

TITANIUM DIOXIDE

2.000

WATER (14,5%) TARGET MOISTURE

**GELATIN** 

COLORANT

TITANIUM DIOXIDE

C.I. NUMBER

77891

E.E.C. NUMBER

E171

### REGULATORY COMPLIANCE

#### **RAW MATERIAL**

GELATIN: Complies with the requirements of current EP and USP/NF editions.

It is of pure bovine origin, complies with revision in force of European Guideline EMEA/410/01. Each supplier has a Ph. Eur certificate of suitability for their product (R1-CEP-2000-027-Rev.02, R1-CEP-2000-029-Rev.05, R1-CEP-2000-045-Rev.03, R1-CEP-2002-110-Rev.00, R1-CEP-2001-211.Rev.01). COLORANTS: Comply with Directive 2009/35, Commission Regulation 231/2012 and where applicable, with the requirements of the EP and USP/NF. PRINTING INKS: Comply with pharmaceutical regulations.

Do not contain preservatives and have not been treated with Ethylene Oxide. They comply with The European Agency for the Evaluation of Medicinal Products (EMEA) guideline CPMP/ICH/283/95 and the European Pharmacopoeia for residual solvents.

DISINTEGRATION: Less than 15 minutes by the Ph, Eur, test and have a rupture time of less than 5 minutes by the USA Federal Specification 285A Acid Solubility Test.

This is to certify that the information above has been approved by QA as conforming to Qualicaps Europe specifications, as described in the current technical manual and applicable regulatory requirements.

**FEDERICO GARCIA** QUALITY ASSURANCE MANAGEMENT

Page

DATE: 31/08/2017 9001-14001